REF UTI-1S1

# Rapid Response® Urinary Tract Infection (UTI) Test Strips For Self Testing

## Instructions for Use


For *in vitro* Diagnostic Use Only.

Unfold leaflet and flip this sheet over for instructions.

## Principle

The Rapid Response® Urinary Tract Infection (UTI) Test Strips are for the detection of leukocytes and nitrite in urine. It is used as an aid in the screening of urinary tract infection (UTI). The test is a firm plastic strip onto which leukocyte and nitrite test pads are attached. If the test is positive, the leukocyte test pad should be beige to dark purple, and the nitrite test pad should be uniform pink to red.

## Composition

| Reagent          | Read Time | Composition  | Description   |
|------------------|-----------|--|---|
| Leukocytes (LEU) | 2 minutes | Derivatized pyrrole amino acid ester; diazonium salt; buffer; non-reactive ingredients | Detects leukocytes as low as 9-15 white blood cells Leu/ $\mu$ L in clinical urine  |
| Nitrite (NIT)    | 1 minutes | P-arsanilic acid; N-(1-naphthyl) ethylenediamine; non-reactive ingredients             | Detects sodium nitrite as low as 0.05-0.1 mg/dL in urine with a low specific gravity and less than 30 mg/dL ascorbic acid |

## Storage and Handling

- Store in a dry place at 2-30°C (36-86°F). Do not freeze. Keep out of direct sunlight.
- The strip should remain in the sealed pouch until use.
- Excessive or prolonged exposure to air before use may cause inaccurate reading.
- Do not transfer the strips to another container.
- Do not use if pouch is torn or damaged.
- Do not touch the test pad(s) of the strip. Discard any discoloured strips that may have deteriorated.

## Warnings and Precautions

Please read all the information in this package insert before performing the test.

- For urine testing only. Do not use for blood testing.
- Keep out of the reach of children.
- The used strip should be discarded according to local regulations after testing.
- Do not use after the expiration date.
- For *in vitro* diagnostic use. Not to be taken internally.

## Limitations

There is the possibility that this test may produce false results. Consult your physician before making any medical decisions.

**Note:** The test may be affected by substances that cause abnormal urine colour such as drugs containing azo dyes (e.g. Pyridium®, Azo Gantrisin®, Azo Gantanol®), nitrofurantoin (Microdantin®, Furadantin®), and riboflavin. The colour development on the test pad may be masked or a colour reaction may be produced that could be interpreted as false results.

1. Leukocytes: The drug tetracycline may cause a false negative result. High protein or elevated glucose in urine may cause test results to be low.
2. Nitrite: Any degree of uniform pink to red colour should be read as a positive result. Pink spots or pink edges should not be read as a positive result. High ascorbic acid may cause test results to be low.

**Note:** You may get a negative result if you have a UTI caused by bacteria that does not change nitrate to nitrite, when urine has not

been held in the bladder for more than 4 hours, when taking antibiotics, or when your diet does not include nitrates.

## Performance Characteristics

The results of the Rapid Response® Urinary Tract Infection (UTI) Test Strips were compared to another commercial test. The study showed 94.4% accuracy for the leukocyte test and 98.4% for the nitrite test.

The performance characteristics of the strips are based on both laboratory and clinical tests. The sensitivity of the test depends upon several factors: the variability of colour perception; the presence or absence of inhibitory factors; and the lighting conditions under which the strip is read.

## FAQ

### What are the common symptoms of UTI?

- The common symptoms of UTI include burning or pain during urination, frequent urination, fever, and back or groin pain.

### Can I have UTI without symptoms?

- Yes. UTI can occur without symptoms in young children, pregnant women and the elderly.

### Do I have to test with first morning urine?

- Although you can test at any time of the day, your first morning urine would have the most bacteria in it if UTI is present. Avoid testing during your menstrual period.

### Do certain foods, alcohol or common medications affect the test?

- In general, most substances will not interfere with the test. Consult your physician if you are taking antibiotics. Some medications (e.g. vitamin B2) can affect the colour of your urine and may cause false results. Large quantities of vitamin C may cause false negative results.

### What if I have questions about how to run the test?

- Please contact support@btnx.com. Medical questions should be discussed with your doctor.

## Glossary of Symbols

Consult instructions for use  
 Test per Kit  
 Store between 2°C to 30°C  
 Use by date  
 Do Not Reuse  
 In vitro diagnostic medical device  
 Lot Number  
 Manufacturer  
 Catalogue number

BTNX, Inc.  
 722 Rosebank Road  
 Pickering, ON L1W 4B2 Canada  
 Technical Support: 1-888-339-9964

# Rapid Response®

## Urinary Tract Infection (UTI) Test Strips For Self Testing

Allow the strip and urine specimen to reach room temperature 15-30°C (59-86°F) prior to testing.

### Setting Up the Test

Gently wash the genital area thoroughly with soap and water then rinse well prior to testing.

### Materials Provided

Unpack the test components from the kit and make sure that all the packaging is intact.

For each test you will need:


Sealed Pouch  
Containing Test Strip  
(with Colour Chart)

Instructions  
(This Document)

Timer or Watch with  
Second Hand  
(Not provided)


Absorbent Material  
(e.g. paper towel)  
(Not provided)


Specimen Collection  
Container  
(Optional)  
(Not provided)

### Instructions


Remove the strip from the sealed pouch and use it as soon as possible.

Hold the end of the strip farthest away from the test pad(s) and begin urinating. After 1 or 2 seconds of passing urine, hold the strip downward with the test pads pointing into the urine stream for 1-2 seconds. Make sure the test pads are completely wet.


OR


**NOTE:** If you prefer, you can urinate into a clean and dry container. Dip the test pad(s) of the strip completely into the urine for 1-2 seconds and remove immediately. Run the edge of the strip against the rim of the urine container to remove excess urine.

After removing the strip from your urine, immediately bring the edge of the strip into contact with an absorbent material (e.g. paper towel) to remove excess urine. Lay the strip with the test pads(s) facing upwards and begin timing.

**Read result at 1 minute for nitrite and at 2 minutes for leukocytes.** Compare the test pad(s) to the colour blocks on the colour chart. Hold the strip close to the colour blocks and carefully match each pad to the colour chart for that test.

**NOTE: Do not read results after 3 minutes.**

### Reading the Results

Read the results by comparing the test pad(s) to the colour blocks on the colour chart. Match the colour of the test pad to the closest colour block on the colour chart.

| Color Chart Reading |         |     | Recommendations  |
|---------------------|---------|-----|--|
| Results             | Reagent |     |  |
|                     | LEU     | NIT |  |
| - or ±              | -       | -   | Results suggest no sign of UTI. When trace (±) LEU occurs or if you have symptoms, test again the following day with your first morning urine using a new strip. If you get another trace (±) LEU or if you still have symptoms, consult your physician. |
| - or ±              | +       | +   | Results suggest a sign of UTI. Consult your physician immediately.   |
| +, ++, or +++       | +       | +   | Results suggest a sign of UTI. Consult your physician immediately.   |
| +, ++, or +++       | -       | -   | Results suggest a sign of UTI. Repeat the test next time you urinate using a new strip. Make sure to wash the genital area first. If LEU is still positive, consult your physician.  |

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BTNX Inc.